

Purification Protocol for His6-Sds3 CCSID (aa 43-234; pMCSG vector) (20 °C overexpression)

Prepare Following Buffer Solutions (sterile filter):

→ Native Buffer:

- *50 mM Tris
- *200 mM NaCl
- *2 mM TCEP
- *pH 8

→ High Salt Buffer:

- *50 mM Tris
- *1 M NaCl
- *2 mM TCEP
- *pH 8

→ 6M Gdn.HCl Buffer:

- *50 mM Tris
- *200 mM NaCl
- *2 mM TCEP
- *6 M Gdn.HCl
- *pH 8

→ Elution Buffer:

- *150 mM EDTA
- *200 mM NaCl
- *2 mM TCEP
- *6 M Gdn.HCl
- *50 mM Tris
- *pH 8

→ Prepare Ni²⁺ column

- Wash resin with 3 column volumes of Gdn.HCl buffer
- Wash resin thoroughly with UV/UF H₂O (10 column volumes)
- Charge resin by adding 5 column volumes 0.1 M Ni²⁺ solution
- Wash resin with UV/UF H₂O
- Wash resin with column buffer (native buffer)
- Wash resin with UV/UF H₂O
- Resin is now ready for use

Protocol:

- 1) Resuspend 500 ml pellet in 30 ml Gdn.HCl buffer (pellet from 1 liter culture = 60 ml buffer)
- 2) Add protease inhibitors (scale accordingly for larger pellets)
 - *300 µl PMSF
 - *300 µl Triton X-100
 - *30 µl leupeptin
 - *30 µl pepstatin
- 3) Sonicate 1 hour – 1 second on, 5 seconds off (preset program number 001)
- 4) Centrifuge 30 minutes @ 12,000 rpm
- 5) Incubate soluble supernatant with prepared Ni²⁺ column
 - *1 hour @ 4 °C*
- 6) Collect flow-through
- 7) Wash with Gdn.HCl buffer (30 ml)^
- 8) Wash with native buffer (30 ml)^
- 9) Wash with high salt buffer (30 ml)^
- 10) Wash with native buffer (30 ml)^
- 11) Elute with elution buffer
 - *Incubate 10 ml with resin for 10-15 minutes @ 4 °C
 - *Collect elution
 - *Repeat 3 times
- 12) Combine elution 1 and 2, filter, inject onto HPLC (program: TORC_45C for column 080408)

Notes:

^ = For a 500 ml or 1 l pellet, use 30 ml washes; increase accordingly for larger pellets

* = Dialysis buttons for gel samples

- Remove 100 µl of each sample and place in dialysis button
- Place in water for minimum 3 hours to remove Gdn.HCl
- Run SDS-PAGE gel